

```

  ⌋ library(knitr) ⌋ opts_hunkset(cache=FALSE,tidy=FALSE,highlight=FALSE)
⌋ opts_hunkset(cache = FALSE, tidy = FALSE, fig.align = "center") ⌋ li-
brary(biogas) ⌋ ⌋ options(width=65) ⌋ NTS: Delete following lines before sub-
mission! ⌋ files ⌋- list.files('../R', full.names = TRUE) ⌋ for(i in files) source(i)
⌋ files ⌋- list.files('../data', full.names = TRUE) ⌋ for(j in files) load(j)

```

Calculating methane and biogas production and production rates using manometric methods

Nanna Løjborg and Sasha D. Hafner (sasha.hafner@eng.au.dk)

September 23, 2019

1 Introduction

Manometric biogas calculation methods are commonly used for calculation of cumulative biogas and methane production and production rates prior to determining methane potential and anaerobic biodegradability of a given substrate. In manometric methods bottle headspace pressure is measured under constant temperature conditions from different techniques such as using pressure manometers or transducers, manometer assisted syringes, or low flow pressure [?]. The measured pressure is equivalent to a certain volume fraction and is converted to biogas and methane (if composition is provided) volume data using manometric calculation methods. Pressure manometers requires a pressure to build up inside the headspace of the bottles. Headspace pressure is typically measured with a portable manometer, which evaluates the pressure relative to atmospheric pressure. Determination of CH₄ volume from volumetric measurements required some data processing. For the purpose of this, an additional function (`cumBgDataPrep()`) was created just to handle data preparation (refer to `cumBgDataPrep.R`). `cumBgDataPrep()` sorts and restructures biogas data from volumetric, manometric, gravimetric, and gas density measurements for further calculation of cumulative biogas and methane production, and production rates.

We developed a function to process manometric measurements only, to simplify functions within the biogas package, making it easier for new biogas package users to work with their data. This document provides a brief description of the manometric biogas calculation function (`cumBgMan()`) for new users. We have assumed that readers are familiar with biogas data collection and R.

2 Overview of the function

`cumBgMan()` is a "high-level" function within the biogas package. The purpose of `cumBgMan()` is to convert pressure data collected in the laboratory to cumulative biogas and CH₄ production and to calculate production rates. The function can handle data from any number of bottles. For simple operations

(e.g. interpolation and standardization of biogas volume) `cumBgMan()` is supported by calls to external "low-level" functions (refer to Section ??, Table ??). The "low-level" functions are straight-forward to use, and details can be found in their individual help files. This document describes the use of `cumBgMan()`.

In general, `cumBg*` functions are compiled of four sections: check arguments, restructuring and sorting data, interpolation if needed, and biogas standardization and calculations. Restructuring and sorting of data and interpolation are handled by the external functions `interp()` and `cumBgDataPrep()`, respectively. From `interp()` gas composition, cumulative biogas production, and other variables can be interpolated to a specified time if required. From `cumBgDataPrep()` "wide" and "long" data structure are restructured to "longcombo" data, which is required for `cumBgMan()` to further calculate cumulative biogas and CH₄ production and production rates. Additionally, data is sorted, headspace is added if provided, and composition data is corrected if it seems to be a percentage. Subsequently, the now restructured and sorted data is standardized in `cumBgVol()` by an external function called `stdVol()`.

Two manometric methods are commonly used to evaluate biogas pressure measurements. Method 1 is based on normalized CH₄ concentrations, whereas method 2 accounts for the actual CH₄ in the bottle headspace. Both methods are available through `cumBgMan()` and results is expected to be independent of method. The examples below describe cumulative biogas calculation on two different datasets. The "long" structured dataset is evaluated using manometric method 1, whereas the "longcombo" dataset is evaluated using manometric method 2, as true methane concentrations are provided. All external functions are within the biogas package

2.1 Function arguments

The arguments for the `cumBgMan` function are:

```
## function (dat, comp = NULL, temp = NULL,
##          interval = TRUE, data.struct = "longcombo",
##          id.name = "id", time.name = "time",
##          dat.name = "pres", comp.name = "xCH4",
##          temp.init = NULL, pres.init = NULL, pres.resid = NULL,
##          rh.resid = NULL, rh.resid.init = 1,
##          headspace = NULL, vol.hs.name = "vol.hs",
##          absolute = TRUE, pres.amb = NULL,
##          cmethod = "removed", imethod = "linear",
##          extrap = FALSE, addt0 = TRUE, showt0 = TRUE,
##          dry = FALSE, empty.name = NULL,
##          std.message = !quiet, check = TRUE,
##          temp.std = getOption("temp.std", as.numeric(NA)),
##          pres.std = getOption("pres.std", as.numeric(NA)),
##          unit.temp = getOption("unit.temp", "C"),
##          unit.pres = getOption("unit.pres", "atm"),
```

```
##          quiet = FALSE)
## NULL
```

Most of the arguments have default values, but to calculate CH₄ production we must provide values for at least `dat` (data frame with pressure measurements), `comp` (data frame with gas composition), and `temp` (biogas temperature)¹. along with the names of a few columns in the input data frames. If `temp` argument is not provided, biogas volumes will not be standardized. Similar, if `comp` argument is not provided, `cumBgMan()` will return calculations on biogas only and no calculations on CH₄.

By default `interval = TRUE` and `data.struct` is set as "longcombo". "wide" and "long" structured data will be restructured to "longcombo" internally by `cumBgDataPrep()`, when specified by the `data.struct` argument. When data are cumulative, the `interval` argument should be set to `FALSE` (refer to Section ??).

Similarly, there is an `id.name` argument for the bottle identification code (ID) column (default = "id"). For `data.struct = "wide"`, there is no ID column. Instead data for each bottle, have individual columns and column names, which are used as ID codes. Here, the name of the column containing the response variables (`dat.name`), is set as the name of the first column with response variables. All following columns are assumed to also have measurement data.

Furthermore, we need to specify the name of the time column containing time data using the `time.name` argument (default = "time"). If separate data frames are used for `dat` and `comp`, the name must match. Time data may be POSIX objects, but then `t0` will not be added to rows by the `cumBgMan` function. In addition the `addt0` argument is used to add row with "time zero" (`time.name = 0`) for each bottle in order to calculate production rates for the first observation (default = `TRUE`). Whereas, `showt0` determines if the "time zero" should be returned in the output (default = `TRUE` if `time.name` is numeric and contains 0 and otherwise `FALSE`).

The `comp.name` argument is used to indicate which column within the `comp` data frame contains CH₄ concentrations (as mole fraction in dry biogas, normalized so the sum of mole fractions of CH₄ and CO₂ sum to unity). Default is `xCH4`. `comp` may also just be a single value instead of a data frame or column. When providing a single value for `comp`, this value is applied to all observation-sand the `comp.name` argument is not needed. The definition of `xCH4` depends on the `cmethod` argument. By default (`cmethod = "removed"`) the function calculates volumes following [?] as the product of standardized volume of biogas removed and normalized CH₄ content. If results should be based on the sum of methane removed and methane remaining in the bottle headspace, `cmethod` should be set to `"total"`. When `cmethod = "total"`, CH₄ concentration is calculated using all components (CH₄, CO₂, N₂, H₂S, etc.) instead of CH₄ and

¹. By default, temperature is in °C and pressure in atm, but these can be changed in the function call with the `temp.unit` and `pres.unit` arguments, or globally with `options`. The same default values apply for temperature and pressure for presentation of biogas and methane, but these can be changed in the function call with the `temp.std` and `pres.std` arguments.

CO₂ only.

If any CH₄ measurements are missing, xCH₄ is interpolated by the external "low-level" function `interp()`. Here, the `imethod` argument can be used to define interpolation method (default = "linear"), which is passed to `interp()`. Similar, an `extrap` argument is passed to `interp()` (default = "FALSE"). The `extrap` argument is used to indicate if composition data (`comp.name`) should be extrapolation (e.g., in the case of missing initial composition values). Extrapolation is constant (the value of the nearest observation is taken). In general, extrapolation and extensive interpolation should be avoided.

Initial headspace temperature and pressure are required to determine initial gas volume in the bottles and can be set using the `temp.init` and `pres.init`, respectively. Default values are NULL. Similarly, post venting headspace pressure is required for manometric calculation methods and can be set using the `pres.resid` argument. If only a single value is provided, this will be used for all observations. Initial and post venting headspace pressure can be absolute or gauge, depending on the `absolute` argument. By default, pressure is absolute (`absolute = TRUE`). When `absolute = FALSE` pressure measurements need to be corrected using the `pres.amb` argument, representing the absolute ambient pressure. `pres.amb` can only be set as a single value and is set as 1 atm or 101.325 kPa by default.

Additionally, a data frame containing headspace volumes is required if `cmethod = "total"` (see description below) and should contain at least a headspace volume column (`vol.hs.name`) and a bottle identification column, with the same column name as in the `dat` and `comp` data frames. The headspace volume column can be set using the `vol.hs.name` argument (default = "vol.hs"). Initial relative humidity of gas in the headspace is set to 1 by default using the `rh.resid.init` argument. If values are provided for relative humidity of gas in the headspace after experiment was initiated, these can be defined by using the `rh.resid` argument.

2.2 Data Structures

Input data may be structured in one of three ways: "long", "wide", or "long-combo". Default is "longcombo", where the composition column is in the `dat` data frame and no composition data frame is required. All three data structure inputs are accepted, but the manometric calculation methods within `cumBgMan()` only process "longcombo" data. "wide" and "long" data are restructured internally by the low-level function `cumBgDataPrep()`. In the following examples all three data structures will be addressed. For more details about each of the three data structure see `biogas_quick_start` vignette.

3 Examples: Calculation of cumulative production of biogas and CH₄ and production rates from manometric methods

Calculation of cumulative biogas and CH₄ production and production rates, typically requires two data frames: Biogas volume measurements and biogas composition (CH₄ fraction)

3.1 "longcombo" data structure

In this example, we will use a "longcombo" example dataset included in the biogas package: `sludgeTwoBiogas` for both biogas volumes and composition and `sludgeTwoSetup` for grouping and headspace volumes. Substrate and inoculum masses provided from `sludgeTwoSetup` are not interesting before calculation of BMP using `summBg()` (refer to Section ??).

These data are from an experiment carried out with seven different substrates all from primary wastewater sludge. The experiment included 24 batch bottles, all with inoculum and one with no substrate (blank):

- Three bottles with inoculum only (Blank50)
- Three bottles with wastewater sludge 25 (WWS25)
- Three bottles with wastewater sludge 25b (WWS25b)
- Three bottles with wastewater sludge 40 (WWS40)
- Three bottles with wastewater sludge 50 (WWS50)
- Three bottles with wastewater sludge 50b (WWS50b)
- Three bottles with wastewater sludge 60 (WWS60)
- Three bottles with wastewater sludge 75 (WWS75)

Data were originally collected by Sergi Astals at the University of Queensland. For more information, please contact Sasha D. Hafner at sasha.hafner@eng.au.dk.

```
> #data("sludgeTwoBiogas")
>
> dim(sludgeTwoBiogas)

[1] 324    8

> head(sludgeTwoBiogas)

  id  time.d pres mass.init mass.final  xCH4  xCO2  xCH4n
1  1  0.0000000    0  187.4168  187.4168  0.0000  0.0000  0.0000000
2  1  0.6270833  368  187.4149  187.3814  0.1010  0.1004  0.5014896
```

```

3 1 1.6027778 286 187.3823 187.3568 0.2253 0.1610 0.5832255
4 1 1.9972222 134 187.3538 187.3447 0.2485 0.1800 0.5799300
5 1 2.7152778 142 187.3406 187.3287 0.3246 0.2062 0.6115298
6 1 3.6180556 168 187.3295 187.3150 0.3674 0.2222 0.6231343

```

```
> summary(sludgeTwoBiogas)
```

```

      id      time.d      pres
Min.   : 1.0   Min.   : 0.000   Min.   : 0.0
1st Qu.: 8.0   1st Qu.: 2.715   1st Qu.: 181.0
Median :12.5   Median : 8.637   Median : 358.0
Mean   :12.0   Mean   :18.113   Mean   : 565.4
3rd Qu.:17.0   3rd Qu.:24.029   3rd Qu.: 753.0
Max.   :21.0   Max.   :83.662   Max.   :3785.0
      mass.init   mass.final   xCH4
Min.   :144.7   Min.   :144.7   Min.   :0.0000
1st Qu.:163.6   1st Qu.:163.6   1st Qu.:0.5479
Median :183.5   Median :183.5   Median :0.6201
Mean   :184.3   Mean   :184.2   Mean   :0.5460
3rd Qu.:200.0   3rd Qu.:199.9   3rd Qu.:0.6418
Max.   :228.4   Max.   :228.4   Max.   :0.7377
      xCO2      xCH4n
Min.   :0.0000   Min.   :0.0000
1st Qu.:0.2599   1st Qu.:0.6446
Median :0.3012   Median :0.6602
Mean   :0.2791   Mean   :0.6220
3rd Qu.:0.3354   3rd Qu.:0.6727
Max.   :0.6344   Max.   :0.7637

```

```
> #data("sludgeTwoSetup")
```

```
>
```

```
> dim(sludgeTwoSetup)
```

```
[1] 18 5
```

```
> head(sludgeTwoSetup)
```

```

      id descrip  vol.hs m.inoc m.sub.vs
1 1 Blank50  79.980  80.02 0.0000000
2 2 Blank50  79.980  80.02 0.0000000
3 3 Blank50  79.950  80.05 0.0000000
4 7  WWS75 117.575  35.02 0.3588463
5 8  WWS75 117.535  35.06 0.3588463
6 9  WWS75 117.490  35.08 0.3600578

```

```
> summary(sludgeTwoSetup)
```

id		descrip	vol.hs	
Min.	: 1.00	Length:18	Min.	: 38.91
1st Qu.:	8.25	Class :character	1st Qu.:	63.19
Median :	12.50	Mode :character	Median :	80.56
Mean :	12.00		Mean :	80.05
3rd Qu.:	16.75		3rd Qu.:	99.43
Max.	:21.00		Max.	:117.58

m.inoc		m.sub.vs	
Min.	: 35.02	Min.	:0.0000
1st Qu.:	50.06	1st Qu.:	0.3591
Median :	72.59	Median :	0.5864
Mean :	68.39	Mean :	0.5605
3rd Qu.:	80.05	3rd Qu.:	0.8130
Max.	:100.08	Max.	:1.0191

Calculating cumulative biogas and CH₄ production and production rates, is the first step in processing data from a BMP trial. Subsequently, BMP can be calculated by the high-level function `summBg()` included in the biogas package. Cumulative biogas and CH₄ production and production rates from pressure data with `sludgeTwoBiogas` and `sludgeTwoSetup` data frames as input, can be calculated from `cumBgMan()`.

To calculate CH₄ production from these "longcombo" data, we must provide values for at least `dat` and `comp`, which is in a combined data frame (we will use `sludgeTwoBiogas`), and `temp` (biogas temperature) along with the names of a few columns in our input data frame.

We can use the default values "longcombo", "id", and "pres" for the `data.struct`, `id.name`, and `dat.name` arguments, respectively. Whereas, the `time.name` and `comp.name` arguments need to be specified as "time.d" and "xCH4n", respectively. Similar, default values can be used for `cmethod` = "removed", evaluating CH₄ concentration based on normalized CH₄ and CO₂ values, for `imethod` = "linear", resulting in internal linear interpolation of xCH4n by calling the `interp` function, and for `extrap` = FALSE.

Initial headspace temperature and pressure and post venting headspace pressure are set as constant values using the `temp.init`, `pres.init`, and `pres.resid` arguments, respectively. Initial and post venting headspace pressure can be absolute or gauge depending on the value of the `absolute` argument. We will set `absolute` = FALSE to account for pressure measurements provided in gauge. In order to calculate absolute pressure from gauge pressure measurements when `absolute` = FALSE, a single absolute ambient pressure value is required. Here we will set `pres.amp` to 1013 mbar. Note the unit of the pressure data in `straw-Pressure`. Default unit is atm, but this can be changed using the `unit.pres` argument. In this example we will set the `unit.pres` = "mbar" to match unit of `dat.name` and `pres.resid` column. Absolute ambient pressure should be provided in the same unit as defined in the `unit.pres` argument.

Finally, headspace volumes are provided from `sludgeTwoSetup`. The data frame containing headspace volumes is defined using the `headspace` argument,

whereas the default value "vol.hs" can be used for the column containing headspace volume data.

```
> cum.prod.lc <- cumBgMan(sludgeTwoBiogas, temp = 30,
+                          time.name = "time.d", comp.name = "xCH4n",
+                          temp.init = 30, pres.init = 0.0, pres.resid = 0,
+                          headspace = sludgeTwoSetup,
+                          pres.amb = 1013, absolute = FALSE,
+                          unit.pres = "mbar")
```

Note the message about standard temperature and pressure—it is important to make sure these values are correct, therefore users are reminded by a message². Also, note warning about pressure unit. This warning is to make sure all pressure measurements are gauge, as absolute is set to FALSE, meaning pressure data are corrected to absolute pressure by the cumBgMan function.

The output looks like this:

```
> head(cum.prod.lc)
```

	id	time.d	pres	mass.init	mass.final	xCH4	xCO2	xCH4n
1	1	0.0000000	0	187.4168	187.4168	0.0000	0.0000	0.0000000
2	1	0.6270833	368	187.4149	187.3814	0.1010	0.1004	0.5014896
3	1	1.6027778	286	187.3823	187.3568	0.2253	0.1610	0.5832255
4	1	1.9972222	134	187.3538	187.3447	0.2485	0.1800	0.5799300
5	1	2.7152778	142	187.3406	187.3287	0.3246	0.2062	0.6115298
6	1	3.6180556	168	187.3295	187.3150	0.3674	0.2222	0.6231343

	vol.hs	temperature	pres.resid	pres.abs	pres.resid.abs
1	79.98	30	0	1013	1013
2	79.98	30	0	1381	1013
3	79.98	30	0	1299	1013
4	79.98	30	0	1147	1013
5	79.98	30	0	1155	1013
6	79.98	30	0	1181	1013

	rh.resid	pres.resid.prev	rh.resid.prev	temp.prev	xCH4n.prev
1	1.0000000	1013	1.0000000	30	0.0000000
2	0.7335264	1013	1.0000000	30	0.0000000
3	0.7798306	1013	0.7335264	30	0.5014896
4	0.8831735	1013	0.7798306	30	0.5832255
5	0.8770563	1013	0.8831735	30	0.5799300
6	0.8577477	1013	0.8770563	30	0.6115298

	vBg	vCH4	cvBg	cvCH4	rvBg	rvCH4
1	0.000000	0.000000	0.000000	0.000000	NA	NA
2	26.168144	13.123052	26.16814	13.12305	41.72993	20.927125
3	19.530015	11.390402	45.69816	24.51345	20.01653	11.674149

²Remember that standard conditions can be set in the function call with `temp.std` and `pres.std`, or globally with `options()`.

```

4  8.861695  5.139163 54.55985 29.65262 22.46627 13.028864
5  9.743608  5.958506 64.30346 35.61112 13.56943  8.298113
6 11.573914  7.212103 75.87738 42.82323 12.82034  7.988791

```

```
> dim(cum.prod.lc)
```

```
[1] 324 24
```

The data frame that is returned has maintained the "longcombo" structure with all the original columns in `sludgeTwoBiogas`, plus additional columns from manometric biogas calculations.

In these columns, `v` stands for (standardized) volume, `cv` (standardized) cumulative volume, `rv` stands for (standardized) volume production rate, and `Bg` and `CH4` for biogas and methane. So `cvBg` contains standardized cumulative biogas production and `cvCH4` contains standardized cumulative CH₄ production.

It is probably easier to understand the data in the output graphically. Here we will use the `ggplot` function from the `ggplot2` package to plot it.

```

> library(ggplot2)
> ggplot(cum.prod.lc, aes(time.d, cvCH4, colour = factor(id))) +
+   geom_point() +
+   geom_line(aes(group = id)) +
+   labs(x = "Time [d]", y = "cvCH4 [mL]", colour = "Bottle id") +
+   theme_bw()

```

3.2 "long" data structure

In this example, we will use a long example data set included in the biogas package: **strawPressure** for headspace pressure in batch bottles, **strawComp** for methane content of biogas, and **strawSetup** for grouping and headspace volumes.

These data are interval-based BMP measurement of headspace pressure in 12 anaerobic bottles with straw as the substrate. The experiment included 12 batch bottles on the following six variables:

- Two bottles with inoculum and straw treated with treatment r3
- Two bottles with inoculum and straw treated with treatment r5
- Two bottles with inoculum and straw treated with treatment r6.5
- Two bottles with inoculum and straw treated with treatment r8
- Two bottles with inoculum and straw treated with treatment r8 no buff
- Two bottles with inoculum and straw treated with treatment r9

Bottles were ca. 600 mL glass serum bottles with butyl rubber septa and screw caps. Pressure was measured using an electronic manometer. Data in **strawMass**, **strawSetup**, and **strawComp** are from the same bottles. More details can be found in the help file for these data.

```
> #data("strawPressure")
>
> dim(strawPressure)
```

```
[1] 72  5
```

```
> head(strawPressure)
```

	bottle		date.time	time	pres	pres.resid
1	1	2016-04-28	10:30:00.00	1.9	153.2034	101.629
2	2	2016-04-28	10:30:00.00	1.9	152.0888	101.629
3	3	2016-04-28	10:30:00.00	1.9	142.0576	101.629
4	4	2016-04-28	10:30:00.00	1.9	142.1590	101.629
5	5	2016-04-28	10:30:00.00	1.9	138.7139	101.629
6	6	2016-04-28	10:30:00.00	1.9	139.5245	101.629

```
> summary(strawPressure)
```

	bottle		date.time	time
Min.	: 1.00	2016-04-28	10:30:00.00:12	Min. : 1.900
1st Qu.:	3.75	2016-04-29	16:00:00.00:12	1st Qu.: 3.120
Median :	6.50	2016-04-30	13:20:00.00:12	Median : 4.965
Mean	: 6.50	2016-05-02	11:00:00.00:12	Mean : 5.987

```

3rd Qu.: 9.25    2016-05-04 14:30:00.00:12    3rd Qu.: 8.060
Max.      :12.00    2016-05-09 10:45:00.00:12    Max.      :12.910
      pres      pres.resid
Min.      :103.9    Min.      :101.0
1st Qu.   :125.3    1st Qu.   :101.6
Median    :138.8    Median    :101.8
Mean      :140.2    Mean      :101.7
3rd Qu.   :155.0    3rd Qu.   :102.0
Max.      :180.7    Max.      :102.1

```

```
> #data("strawComp")
```

```
>
```

```
> dim(strawComp)
```

```
[1] 63  4
```

```
> head(strawComp)
```

	bottle		date.time	time	xCH4
1	1	2016-04-28	10:30:00.00	1.9	0.4189
2	2	2016-04-28	10:30:00.00	1.9	0.4219
3	3	2016-04-28	10:30:00.00	1.9	0.3395
4	4	2016-04-28	10:30:00.00	1.9	0.3015
5	5	2016-04-28	10:30:00.00	1.9	0.2769
6	6	2016-04-28	10:30:00.00	1.9	0.2850

```
> summary(strawComp)
```

	bottle		date.time		time
Min.	: 1.000	2016-04-28	10:30:00.00:12	Min.	: 1.900
1st Qu.	: 3.000	2016-05-04	14:30:00.00:11	1st Qu.	: 3.120
Median	: 6.000	2016-04-30	13:20:00.00:10	Median	: 4.010
Mean	: 6.175	2016-05-09	10:45:00.00: 9	Mean	: 5.724
3rd Qu.	: 9.000	2016-04-29	15:35:00.00: 1	3rd Qu.	: 8.060
Max.	:12.000	2016-04-29	15:40:00.00: 1	Max.	:12.910
		(Other)			:19

```

      xCH4
Min.      :0.01187
1st Qu.   :0.30020
Median    :0.43490
Mean      :0.39450
3rd Qu.   :0.48925
Max.      :0.88720

```

```
> #data("strawSetup")
```

```
>
```

```
> dim(strawSetup)
```

```
[1] 12 6
```

```
> head(strawSetup)
```

	bottle	treatment	start	sub.mass	inoc.mass
1	1	r3	2016-04-26 13:00:00.00	5.27	45.36
2	2	r3	2016-04-26 13:00:00.00	5.27	43.88
3	3	r5	2016-04-26 13:00:00.00	5.28	26.27
4	4	r5	2016-04-26 13:00:00.00	5.30	26.95
5	5	r6.5	2016-04-26 13:00:00.00	5.31	20.37
6	6	r6.5	2016-04-26 13:00:00.00	5.29	20.61

	headspace
1	491.3
2	493.3
3	502.9
4	502.4
5	502.4
6	505.8

```
> summary(strawSetup)
```

	bottle	treatment	start
Min.	: 1.00	r3	:2 2016-04-26 13:00:00.00:12
1st Qu.:	3.75	r5	:2
Median :	6.50	r6.5	:2
Mean :	6.50	r8	:2
3rd Qu.:	9.25	r8 no buff:	2
Max.	:12.00	r9	:2

	sub.mass	inoc.mass	headspace
Min.	:5.240	Min. :14.57	Min. :491.3
1st Qu.:	5.270	1st Qu.:16.48	1st Qu.:502.4
Median :	5.280	Median :18.76	Median :502.9
Mean :	5.284	Mean :23.31	Mean :503.1
3rd Qu.:	5.300	3rd Qu.:26.44	3rd Qu.:507.5
Max.	:5.340	Max. :45.36	Max. :509.5

As with the "longcombo" data, cumulative production of CH₄ is needed in order to calculate BMP. Again, we can calculate these with the `cumBgMan` function, using `strawPressure`, `strawComp`, and `strawSetup` data frames as input.

To calculate CH₄ production from these "long" structured data, we must provide values for at least `dat`, `comp`, and `temp` along with the names of a few columns in our input data frame. The `dat` and `comp` arguments are set as the data frames: `strawPressure` and `strawComp`, whereas `temp` is set as single values of 31 °C. For `data.struct != "longcombo"` the data structure needs to be specified. Here we set `data.struct = "long"` which is internally passed to `cumBgDataPrep()` and restructured to "longcombo" structure prior to entering

the manometric calculations Furthermore, we need to specify the name of the id column in `strawPressure` as `bottle` using the `id.name` argument.

We can use the default values `"time"`, `"pres"`, and `"xCH4"` for the `time.name`, `dat.name` and `comp.name` arguments, respectively. The id and time columns create a link between the two data frames. If observations are not from the same time this is solved by the external `interp` function by interpolation. Similar, default values can be used for `cmethod = "removed"`, evaluating CH₄ concentration based on normalized CH₄ and CO₂ values and for `imethod = "linear"`, resulting in internal linear interpolation of xCH₄ by calling the `interp` function. Additionally, the response variables are interval data only and hence, we can use the default `interval = TRUE`.

In addition to interpolation for later observations, an extrapolation argument (`extrap`) can be provided if required. We do not have initial biogas composition (compare the heads of `strawPressure` and `strawComp`) so extrapolation is required to calculate initial volumes and rates. Therefore, we need to set `extrap = TRUE`.

Initial headspace temperature and pressure are set as constant values using the `temp.init` and `pres.init` arguments, respectively. Whereas, headspace pressure after venting is provided from the `strawPressure` data frame and can be defined by assigning the column name to the `pres.resid` argument. Similar to the "longcombo" data example, we will set `absolute = FALSE` to account for pressure measurements provided in gauge. In order to calculate absolute pressure from gauge pressure measurements when `absolute = FALSE`, a single absolute ambient pressure value is required. To correct gauge pressure data to absolute pressure, we will set `pres.amb` to 101.3 kPa. Again, note the unit of the pressure data in `strawPressure`. Pressure unit is set as atm by default, but this can be changed using the `unit.pres` argument. In this example we will set the `unit.pres = "kpa"` to match unit of `dat.name` and `pres.resid` column.

Finally, headspace volumes are provided from `strawSetup`. The data frame containing these volumes is defined using the `headspace` argument, whereas the column is set as `vol.hs.name = "headspace"`.

```
> cum.prod.1 <- cumBgMan(strawPressure, comp = strawComp, temp = 31,
+                         data.struct = "long",
+                         id.name = "bottle",
+                         temp.init = 21.55, pres.resid = "pres.resid", pres.init = 0.0,
+                         headspace = strawSetup, vol.hs.name = "headspace",
+                         pres.amb = 101.3, absolute = FALSE,
+                         extrap = TRUE,
+                         unit.pres = "kPa")
> head(cum.prod.1)
```

	bottle		date.time	time	pres	pres.resid
1	1		<NA>	0.00	NA	NA
2	1	2016-04-28	10:30:00.00	1.90	153.2034	101.6290
3	1	2016-04-29	16:00:00.00	3.12	171.9485	101.8316

```

4      1 2016-04-30 13:20:00.00 4.01 140.7404 101.8316
5      1 2016-05-02 11:00:00.00 5.92 168.7061 102.1356
6      1 2016-05-04 14:30:00.00 8.06 168.4022 102.0343
      xCH4 headspace temperature pres.abs pres.resid.abs
1      NA      NA      NA      NA      NA
2 0.4189000      491.3      31 254.5034      202.9290
3 0.4363556      491.3      31 273.2485      203.1316
4 0.4413000      491.3      31 242.0404      203.1316
5 0.5014287      491.3      31 270.0061      203.4356
6 0.5152000      491.3      31 269.7022      203.3343
      rh.resid pres.resid.prev rh.resid.prev temp.prev xCH4.prev
1      NA      NA      NA      NA      NA
2 0.7973527      101.3000      1.0000000      21.55 0.0000000
3 0.7433951      202.9290      0.7973527      31.00 0.4189000
4 0.8392467      203.1316      0.7433951      31.00 0.4363556
5 0.7534481      203.1316      0.8392467      31.00 0.4413000
6 0.7539216      203.4356      0.7534481      31.00 0.5014287
      vBg      vCH4      cvBg      cvCH4      rvBg      rvCH4
1 0.0000 0.00000 0.0000 0.0000      NA      NA
2 644.8100 270.11092 644.8100 270.1109 339.3737 142.16364
3 302.1496 131.84467 946.9597 401.9556 247.6636 108.06940
4 164.3333 72.52027 1111.2929 474.4759 184.6441 81.48345
5 287.9775 144.40021 1399.2705 618.8761 150.7736 75.60220
6 283.6379 146.13027 1682.9084 765.0063 132.5411 68.28517

```

Note the message about standard temperature and pressure—it is important to make sure these values are correct, therefore users are reminded by a message³. Also, note warning concerning pressure ranges. The pressure ranges from 101.3 to 203.4 kPa is reasonable and warning can be ignored. We are aware of this issue and it will be addressed (refer to Section ??).

The output looks like this:

```
> head(cum.prod.1)
```

```

      bottle      date.time time      pres pres.resid
1      1      <NA> 0.00      NA      NA
2      1 2016-04-28 10:30:00.00 1.90 153.2034 101.6290
3      1 2016-04-29 16:00:00.00 3.12 171.9485 101.8316
4      1 2016-04-30 13:20:00.00 4.01 140.7404 101.8316
5      1 2016-05-02 11:00:00.00 5.92 168.7061 102.1356
6      1 2016-05-04 14:30:00.00 8.06 168.4022 102.0343
      xCH4 headspace temperature pres.abs pres.resid.abs
1      NA      NA      NA      NA      NA
2 0.4189000      491.3      31 254.5034      202.9290

```

³Remember that standard conditions can be set in the function call with `temp.std` and `pres.std`, or globally with `options()`.

```

3 0.4363556      491.3          31 273.2485      203.1316
4 0.4413000      491.3          31 242.0404      203.1316
5 0.5014287      491.3          31 270.0061      203.4356
6 0.5152000      491.3          31 269.7022      203.3343
  rh.resid pres.resid.prev rh.resid.prev temp.prev xCH4.prev
1      NA              NA              NA      NA      NA
2 0.7973527      101.3000      1.0000000      21.55 0.0000000
3 0.7433951      202.9290      0.7973527      31.00 0.4189000
4 0.8392467      203.1316      0.7433951      31.00 0.4363556
5 0.7534481      203.1316      0.8392467      31.00 0.4413000
6 0.7539216      203.4356      0.7534481      31.00 0.5014287
      vBg      vCH4      cvBg      cvCH4      rvBg      rvCH4
1  0.0000  0.00000  0.0000  0.0000      NA      NA
2 644.8100 270.11092 644.8100 270.1109 339.3737 142.16364
3 302.1496 131.84467 946.9597 401.9556 247.6636 108.06940
4 164.3333  72.52027 1111.2929 474.4759 184.6441  81.48345
5 287.9775 144.40021 1399.2705 618.8761 150.7736  75.60220
6 283.6379 146.13027 1682.9084 765.0063 132.5411  68.28517

> dim(cum.prod.1)

[1] 84 21

```

The data frame that is returned has been restructured to "longcombo" structure and contains all the original columns in **strawPressure**, plus additional columns from volumetric biogas calculations (refer to Section ??)

As with the "longcombo" data example, the results may be easier to understand graphically using the **ggplot** function from the **ggplot2** package.

```

> ggplot(cum.prod.1, aes(time, cvCH4, colour = factor(bottle))) +
+   geom_point() +
+   geom_line(aes(group = bottle)) +
+   labs(x = "Time [d]", y = "cvCH4 [mL]", colour = "Bottle id") +
+   theme_bw()

```


3.3 Additional examples

For examples using "wide" structured input data frame and varying settings of other arguments (e.g. `addt0` and `showt0`), refer to vignette for volumetric `cumBg` function (`cumBgVol_function.Rnw`).

4 Continuing with the `cumBgMan` function

The `cumBgMan` function is one of several `cumBg*` functions within the `bio-gas` package. Results from the `cumBg*` functions can be used directly in the `summBg` function from the `biogas` package to calculate BMP for the data in question. Though, this operation often requires additional setup information (e.g. inoculum and substrate mass), which is most commonly provided in an external data frame (`setup`). More details can be found in the associated help file

In the current version of `cumBgMan()`, an issue is commonly observed with the warning regarding pressure range during cumulative biogas calculations. The warning is sometimes returned with the results even though range is reasonable. This has been raised as an issue and should be addressed in a future scope.

5 Function internals

Table 1: Operations done with the low-level functions in `cumBgMan()`. All functions are vectorized. See help files for more details.

Operation	Function
Standardize gas volume	<code>stdVol()</code>
Interpolate composition etc.	<code>interp()</code>
Structurize and sort data	<code>cumBgDataPrep()</code>

References

- [1] J. Filer, H. H. Ding and S. Chang Biochemical Methane Potential (BMP) Assay Method for Anaerobic Digestion Research. *Water*, 11, 921, 2019.
- [2] B.K. Richards, R.J. Cummings, T.E. White, and W.J. Jewell. Methods for kinetic-analysis of methane fermentation in high solids biomass digesters. *Biomass & Bioenergy*, 1(2):65–73, 1991.